**Introduction**  
Thank you for purchasing the Naipo Foot Massager with Rolling Massage and Air Compression. Please read this user manual carefully and keep it in a safe place for future reference. If you need any assistance, please contact our support team with your product model number and Amazon order number.

**Included**  
1 x Foot massager  
1 (pair) x Insoles  
1 x Adapter  
1 x User manual

**Product Specifications**

Model	MGF-F15
Name	Naipo Foot Massager with Rolling Massage and Air Compression
Input	12V 3A
Power consumption	36W
Timer	10 minutes / 15 minutes / 20 minutes

**Important Safeguards**

1. Please do not use the product in a humid environment to avoid electric shocks and other possible accidents.
2. Please do not use the product in a dusty environment, otherwise malfunction might occur.
3. Please do not expose the product to the sunlight.
4. Please do not sit on the product.
5. People with health problems should consult a doctor before using the product.
6. Supervision from adults is necessary if users are children or people with reduced ability.
7. Pregnant women should not use the product.
8. Please do not use the product if you are to sleep soon.
9. Please stop using the product if severe pain or discomfort is felt during the use.
10. Please do not use the product if it works abnormally.
11. It's recommended that massage time not exceed 20 minutes at a time.

**Q & A**

Problem	Solution
The product does not work / the product stops working.	Please check whether the adapter is connected tightly to the power supply. Please check if the adapter is connected to the product properly. The product shuts off automatically after reaching the preset time. Please restart the product. Please contact the customer service.
The heating effect is not that obvious.	Please wait for 5 minutes after the heating function has been activated.

**Product Diagram**

**1. Power Button**  
1) After connected to the power supply, the product will shift to the standby mode with an audible signal, and the logo will light up [blue].  
2) Press the button to turn on the product, and all functions will be activated by default.

**2. Heating Button**  
When the product is turned on, the heating function will be activated by default (indicator: blue). Press the button to shut off the function with the indicator off.

**3. Intensity Button**  
1) Weak: When the product is turned on, the intensity level will be weak by default (indicator: blue).  
2) Medium (indicator: green)  
3) High (indicator: red)

**4. Mode Button**  
1) Mode 1: Rolling massage + air compression (massage time: rolling massage > air massage / indicator: blue)  
2) Mode 2: Rolling massage + air compression (massage time: air massage > rolling massage / indicator: red)  
3) Mode 3: Air massage (indicator flickers between blue and red alternately)

**5. Timer Button**  
1) 10 minutes: When the product is turned on, the working time of the product will be 10 minutes by default (indicator: blue).  
2) 15 minutes: Press the button, and the working time of the product will be 15 minutes (indicator: red).  
3) 20 minutes: Press the button again, and the working time of the product will be 20 minutes (indicator: purple).

**Product Maintenance**

1. Please disconnect the product from the power supply before cleaning.
2. Please wipe the surface of the product with mildly damp cloth or sponge for cleaning.
3. Please do not clean with corrosive detergent.
4. The foot sleeves could be detached for cleaning.
5. Please store the product in a cool and dry place.

**Warranty & Customer Support**  
For questions, assistance, support, or warranty claims, please contact us via the email or the address below.  
Web: www.naipocare.com  
E-mail: support@naipocare.com  
Address: Huanan City Electronic Trading Building P07 Room 120, Longgang District, Shenzhen, China.

\*Note: Products purchased directly from Naipo benefit from our after-sales service. If the purchases are made from other selling teams, we recommend you contact them for warranty issues.

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